

## Anti-CD300E antibody

<b>Cat. No.</b>	ml224461
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-CD300E rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Fusion protein of human CD300E
<b>Reactivity</b>	Human, Mouse
<b>Content</b>	0.4 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	CD300E
<b>Full name</b>	CD300e molecule
<b>Synonyms</b>	CLM2; CLM-2; IREM2; PIgR2; IREM-2; PIgR-2; CD300LE; CMRF35-A5
<b>Swissprot</b>	Q496F6

### Target Background

This gene encodes a member of the CD300 glycoprotein family of cell surface proteins expressed on myeloid cells. The protein interacts with the TYRO protein tyrosine kinase-binding protein and is thought to act as an activating receptor.

订购热线: 4008-898-798

**Applications**

**Western blotting**

Predicted band size: 23 kDa

Positive control: Mouse lung tissue lysate

Recommended dilution: 500-2000

Gel: 12% SDS-PAGE


Lysate: 40  $\mu$ g

Lane: Mouse lung tissue lysate

Primary antibody: ml224461 (CD300E Antibody) at dilution 1/550

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



**ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)