

Anti-C11orf49 antibody

Cat. No.	ml224816
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-C11orf49 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human C11orf49
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	C11orf49
Full name	chromosome 11 open reading frame 49
Synonyms	
Swissprot	Q9H6J7

Target Background

The function of this protein remains unknown.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 38 kDa

Positive control: Human stomach tissue lysates

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Human stomach tissue lysates

Primary antibody: ml224816(C11orf49 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn