

Anti-FIS1 antibody

Cat. No.	ml222377
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-FIS1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human FIS1
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	FIS1
Full name	fission 1 (mitochondrial outer membrane) homolog (S. cerevisiae)
Synonyms	TTC11; CGI-135
Swissprot	Q9Y3D6

Target Background

The balance between fission and fusion regulates the morphology of mitochondria. TTC11 is a component of a mitochondrial complex that promotes mitochondrial fission.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 17 kDa

Positive control: Raji and 293T cells, mouse spleen tissue and A431 cells

Recommended dilution: 500-2000

Gel: 10% SDS-PAGE


Lysate: 40 µg

Lane 1-4: Raji cells, 293T cells, mouse spleen tissue, A431 cells

Primary antibody: ml222377(FIS1 Antibody) at dilution 1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn