

Anti-FMNL1 antibody

Cat. No.	ml222383
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-FMNL1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human FMNL1
Reactivity	Human, Mouse
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	FMNL1
Full name	formin-like 1
Synonyms	FMNL; FHOD4; KW-13; C17orf1; C17orf1B
Swissprot	O95466

Target Background

This gene encodes a formin-related protein. Formin-related proteins have been implicated in morphogenesis, cytokinesis, and cell polarity. An alternative splice variant has been described but its full length sequence has not been determined.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 122 kDa

Positive control: A172 cells

Recommended dilution: 500-2000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane: A172 cells

Primary antibody: ml222383 (FMNL1 Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 40 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn