

订购热线: 4008-898-798

Anti-NEDD9 antibody

Cat. No. ml222470

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-NEDD9 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human NEDD9

ReactivityHumanContent0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol NEDD9

Full name neural precursor cell expressed, developmentally down-regulated 9

Synonyms CAS2; CASL; HEF1; CAS-L; CASS2

Swissprot Q14511

Target Background

The protein encoded by this gene is a member of the CRK-associated substrates family. Members of this family are adhesion docking molecules that mediate protein-protein interactions for signal transduction pathways. This protein is a focal adhesion protein that acts as a scaffold to regulate signaling complexes important in cell attachment, migration and invasion as well as apoptosis and the cell cycle. This protein has also been reported to have a role in cancer metastasis. Alternative splicing results in multiple transcript variants.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:93 kDa Positive control:Raji cells

Recommended dilution: 200-1000


Gel: 6%SDS-PAGE Lysate: 40 µg Lane: Raji cells

Primary antibody: ml222470(NEDD9 Antibody) at dilution

1/650

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn