

订购热线: 4008-898-798

Anti-GNAI3 antibody

Cat. No. ml122401

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-GNAI3 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human GNAI3

Reactivity Human, Mouse, Rat

Content0.3 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol GNAI3

Full name guanine nucleotide binding protein (G protein), alpha inhibiting activity polypeptide 3

Synonyms 87U6; ARCND1

Swissprot P08754

Target Background

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling pathways. G proteins are composed of 3 units: alpha, beta and gamma. This gene encodes an alpha subunit and belongs to the G-alpha family. Mutation in this gene, resulting in a gly40-to-arg substitution, is associated with auriculocondylar syndrome, and shown to affect downstream targets in the G protein-coupled endothelin receptor pathway.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human lung cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml122401(GNAI3 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn