

Anti-CDK18 antibody

Cat. No.	ml225074
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CDK18 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human CDK18
Reactivity	Human, Mouse, Rat
Content	0.66 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CDK18
Full name	cyclin dependent kinase 18
Synonyms	PCTK3; PCTAIRE; PCTAIRE3
Swissprot	Q07002

Target Background

May play a role in signal transduction cascades in terminally differentiated cells.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 54 kDa

Positive control: HL-60 cell lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: HL-60 cell lysate

Primary antibody: ml225074(CDK18 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn