

Anti-MRM1 antibody

Cat. No.	ml225255
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MRM1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human MRM1
Reactivity	Human
Content	1.02 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MRM1
Full name	mitochondrial rRNA methyltransferase 1
Synonyms	
Swissprot	Q6IN84

Target Background

Probably methylates the ribose of guanosine G-2270 in the peptidyl transferase center of the mitochondrial large ribosomal RNA (21S).

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 39 kDa

Positive control: K562 and HEPG2 cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-2: K562 and HEPG2 cell lysates

Primary antibody: ml225255(MRM1 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn