

Anti-TIMM22 antibody

Cat. No.	ml225301
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TIMM22 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human TIMM22
Reactivity	Human, Mouse, Rat
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TIMM22
Full name	translocase of inner mitochondrial membrane 22
Synonyms	TEX4; TIM22
Swissprot	Q9Y584

Target Background

Multipass transmembrane proteins are brought into mitochondria and inserted into the mitochondrial inner membrane by way of the TIM22 complex. This complex has six subunits and is a twin-pore translocase. The protein encoded by this gene is a subunit of TIM22 and represents the voltage-activated and signal-gated channel.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 20 kDa

Positive control: HEPG2 cell lysate

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane: HEPG2 cell lysate

Primary antibody: ml225301 (TIMM22 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn