

Anti-DTYMK antibody

Cat. No.	ml225320
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-DTYMK rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human DTYMK
Reactivity	Human, Mouse
Content	1.02 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	DTYMK
Full name	deoxythymidylate kinase
Synonyms	CDC8; TMPK; TYMK; PP3731
Swissprot	P23919

Target Background

Catalyzes the conversion of dTMP to dTDP.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 24 kDa

Positive control: Human fetal liver tissue, HepG2, K562 and Raji cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-4: Human fetal liver tissue, HepG2, K562 and Raji cell lysates

Primary antibody: ml225320(DTYMK Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn