

Anti-NIPSNAP3A antibody

Cat. No.	ml225357
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NIPSNAP3A rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human NIPSNAP3A
Reactivity	Human, Mouse
Content	1.8 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NIPSNAP3A
Full name	nipsnap homolog 3A
Synonyms	TASSC; HSPC299; NIPSNAP4
Swissprot	Q9UFN0

Target Background

NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 28 kDa

Positive control: Mouse brain tissue lysate

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Mouse brain tissue lysate

Primary antibody: ml225357(NIPSNAP3A Antibody) at
dilution 1/550

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 3 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn