

Anti-SAMSN1 antibody

Cat. No.	ml225815
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SAMSN1 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human SAMSN1
Reactivity	Human, Mouse
Content	1.38 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SAMSN1
Full name	SAM domain, SH3 domain and nuclear localization signals 1

Synonyms SLy2; HACs1; NASH1; SASH2; SH3D6B

Swissprot Q9NSI8

Target Background

SAMSN1 is a member of a novel gene family of putative adaptors and scaffold proteins containing SH3 and SAM (sterile alpha motif) domains (Claudio et al., 2001 [PubMed 11536050]).


Applications

Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human esophagus cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml225815(SAMS1 Antibody) at dilution 1/80, on the right is treated with fusion protein.
(Original magnification: ×200)

Western blotting

Predicted band size: 42 kDa

Positive control: Mouse liver tissue lysate

Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane: Mouse liver tissue lysate

Primary antibody: ml225815(SAMSN1 Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 1 second



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn