

Anti-AKAP17A antibody

Cat. No.	ml225885
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-AKAP17A rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human AKAP17A
Reactivity	Human
Content	0.84 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	AKAP17A
Full name	A-kinase anchoring protein 17A
Synonyms	XE7; 721P; XE7Y; CCDC133; CXYorf3; PRKA17A; SFRS17A; AKAP-17A; DXYS155E
Swissprot	Q02040

Target Background

This locus encodes a protein kinase A anchoring protein. The encoded protein is part of the spliceosome complex and is involved in the regulation of alternate splicing in some mRNA precursors. Alternatively spliced transcript variants have been identified for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 81 kDa

Positive control: K562, Jurkat, Hela and 293T cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-4: K562, Jurkat, Hela and 293T cell lysates

Primary antibody: ml225885(AKAP17A Antibody)
at dilution 1/650

Secondary antibody: Goat anti rabbit IgG at 1/5000
dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn