

Anti-DPYSL4 antibody

Cat. No.	ml260205
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-DPYSL4 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human DPYSL4
Reactivity	Human, Mouse, Rat
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	DPYSL4
Full name	dihydropyrimidinase-like 4
Synonyms	CRMP3, DRP-4, ULIP4
Swissprot	O14531

Target Background

Collapsin response mediator protein 3 is an enzyme that in humans is encoded by the CRMP3 gene. It is necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. And plays a role in axon guidance, neuronal growth cone collapse and cell migration.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 62 kDa

Positive control: Human fetal brain tissue

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 30 µg

Lane: Human fetal brain tissue

Primary antibody: ml260205(DPYSL4 Antibody) at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 500-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn