

订购热线: 4008-898-798

Anti-EIF4G2 antibody

Cat. No. ml260254

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-EIF4G2 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human EIF4G2

Reactivity Human, Mouse
Content 0.4 mg/ml
Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol EIF4G2

Full name Eukaryotic translation initiation factor 4 gamma 2

Synonyms P97, AAG1, DAP5, NAT1

Swissprot P78344

Target Background

Translation initiation is mediated by specific recognition of the cap structure by eukaryotic translation initiation factor 4F (eIF4F), which is a cap binding protein complex that consists of three subunits: eIF4A, eIF4E and eIF4G. The protein encoded by this gene shares similarity with the C-terminal region of eIF4G that contains the binding sites for eIF4A and eIF3; eIF4G, in addition, contains a binding site for eIF4E at the N-terminus. Unlike eIF4G, which supports cap-dependent and independent translation, this gene product functions as a general repressor of translation by forming translationally inactive complexes. In vitro and in vivo studies indicate that translation of this mRNA initiates exclusively at a non-AUG (GUG) codon. Alternatively spliced transcript variants encoding different isoforms of this gene have been described.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:103 kDa

Positive control:Hela, HepG2 and NIH/3T3 cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE Lysate: 40 μg


Lane 1-3: Hela cells, HepG2 cells, NIH/3T3 cells

Primary antibody: ml260254(EIF4G2 Antibody) at dilution

1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn