

订购热线: 4008-898-798

Anti-SFRP1 antibody

Cat. No. ml260294

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-SFRP1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human SFRP1

ReactivityHuman, MouseContent0.72 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol SFRP1

Full namesecreted frizzled-related protein 1SynonymsFRP; FRP1; FrzA; FRP-1; SARP2

Swissprot Q8N474

Target Background

This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. Members of this family act as soluble modulators of Wnt signaling; epigenetic silencing of SFRP genes leads to deregulated activation of the Wnt-pathway which is associated with cancer. This gene may also be involved in determining the polarity of photoreceptor cells in the retina.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:35 kDa
Positive control:293T cell lysate
Recommended dilution: 200-1000

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane: 293T cell lysate

Primary antibody: ml260294(SFRP1 Antibody) at dilution

1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn