

订购热线: 4008-898-798

Anti-PTGER4 antibody

Cat. No. ml260305

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-PTGER4 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human PTGER4

ReactivityHumanContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol PTGER4

Full name prostaglandin E receptor 4

Synonyms EP4; EP4R Swissprot P35408

Target Background

The protein encoded by this gene is a member of the G-protein coupled receptor family. This protein is one of four receptors identified for prostaglandin E2 (PGE2). This receptor can activate T-cell factor signaling. It has been shown to mediate PGE2 induced expression of early growth response 1 (EGR1), regulate the level and stability of cyclooxygenase-2 mRNA, and lead to the phosphorylation of glycogen synthase kinase-3. Knockout studies in mice suggest that this receptor may be involved in the neonatal adaptation of circulatory system, osteoporosis, as well as initiation of skin immune responses.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:53 kDa Positive control:Raji cells

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg Lane: Raji cells

Primary antibody: ml260305(PTGER4 Antibody) at dilution

1/450

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn