

Anti-FOXD1 antibody

Cat. No.	ml260349
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-FOXD1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human FOXD1
Reactivity	Human, Mouse
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	FOXD1
Full name	Forkhead box protein D1
Synonyms	FKHL8, FREAC4
Swissprot	Q16676

Target Background

This intronless gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in tumor formation.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 46 kDa

Positive control: A549 cells

Recommended dilution: 500-2000

Gel: 8% SDS-PAGE


Lysate: 40 µg

Lane: A549 cells

Primary antibody: ml260349 (FOX D1 Antibody) at dilution 1/670

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn