

Anti-CMTM5 antibody

Cat. No.	ml260395
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CMTM5 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CMTM5
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CMTM5
Full name	CKLF-like MARVEL transmembrane domain containing 5
Synonyms	CKLFSF5
Swissprot	Q96DZ9

Target Background

This gene belongs to the chemokine-like factor gene superfamily. This family of genes encodes multi-pass membrane proteins that are similar to both the chemokine and the transmembrane 4 superfamilies of signaling molecules. Alternate transcriptional splice variants of this gene, encoding different isoforms, have been characterized. Highly expressed in the brain.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 13 kDa

Positive control: Mouse heart tissue

Recommended dilution: 500-2000

Gel: 12+15% SDS-PAGE


Lysate: 40 µg

Lane: Mouse heart tissue

Primary antibody: ml260395 (CMTM5 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn