

Anti-MC5R antibody

Cat. No.	ml260592
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MC5R rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human MC5R
Reactivity	Human, Mouse
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MC5R
Full name	melanocortin 5 receptor
Synonyms	MC2
Swissprot	P33032

Target Background

This gene encodes a member of the seven-pass transmembrane G protein-coupled melanocortin receptor protein family that stimulate cAMP signal transduction. The encoded protein is a receptor for melanocyte-stimulating hormone and adrenocorticotrophic hormone and is suggested to play a role in sebum generation.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 37 kDa

Positive control: Mouse muscle tissue

Recommended dilution: 500-2000

Gel: 12% SDS-PAGE


Lysate: 30 µg

Lane: Mouse muscle tissue

Primary antibody: ml260592(MC5R Antibody) at dilution 1/1200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 40 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn