

Anti-NCAN antibody

Cat. No.	ml260667
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NCAN rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human NCAN
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NCAN
Full name	neurocan
Synonyms	CSPG3
Swissprot	O14594

Target Background

Neurocan is a chondroitin sulfate proteoglycan thought to be involved in the modulation of cell adhesion and migration. May modulate neuronal adhesion and neurite growth during development by binding to neural cell adhesion molecules (NG-CAM and N-CAM). Chondroitin sulfate proteoglycan; binds to hyaluronic acid.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 143 kDa

Positive control: Mouse brain tissue

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Mouse brain tissue

Primary antibody: ml260667(NCAN Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn