

Anti-ROBO4 antibody

Cat. No.	ml123052
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ROBO4 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human ROBO4
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ROBO4
Full name	roundabout, axon guidance receptor, homolog 4 (Drosophila)
Synonyms	MRB; ECSM4
Swissprot	Q8WZ75

Target Background

Robo4 (roundabout homolog 4 or magic roundabout) is a receptor for Slit proteins which implicate in angiogenesis and vascular patterning especially through facilitating inhibition of primary endothelial cell migration. It is a fairly recently identified member of Robo family and is shorter than other Robo members at both intracellular as well as extracellular domains. Robo4 plays important roles in several aspects of vascular development such as the guidance of endothelial cell migration, cell cycle regulation, inhibition of pathologic angiogenesis and endothelial hyperpermeability etc. Employing mouse and zebrafish models, Slit2/Robo4 interactions have been shown to maintain integrity of vascular network and its barrier function by inhibiting cytokine-mediated vasculogenesis with enhanced permeability, and that Slit2-Robo4-paxillin-GIT1 network inhibits neovascularization and vascular leakage. Moreover, single nucleotide polymorphism of Robo4 has been documented to be associated with autism suggesting its role in development and function of human brain.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human brain

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml123052(ROBO4 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn