

Anti-KDM6B antibody

Cat. No.	ml260900
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-KDM6B rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human KDM6B
Reactivity	Human, Mouse
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	KDM6B
Full name	lysine (K)-specific demethylase 6B
Synonyms	JMJD3
Swissprot	O15054

Target Background

KDM6B (Lysine demethylase 6B) is a 1,679 amino acid nuclear protein that contains one JmjC domain and belongs to the highly conserved JmjC domain-containing protein family. Functioning as a histone demethylase, JMJD3 uses iron and ascorbate as cofactors to demethylate dimethylated and trimethylated Lys-27 residues of Histone H3, thereby playing an important role in the modification of the histone code.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:177 kDa

Positive control:A172 cells

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE

Lysate: 40 µg

Lane: A172 cells

Primary antibody: ml260900(KDM6B Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes



ELISA

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn