

订购热线: 4008-898-798

Anti-NOX5 antibody

Cat. No. ml222775

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-NOX5 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human NOX5

ReactivityHumanContent0.1 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol NOX5

Full name NADPH oxidase, EF-hand calcium binding domain 5

Synonyms

Swissprot Q96PH1

Target Background

This gene is predominantly expressed in the testis and lymphocyte-rich areas of spleen and lymph nodes. It encodes a calcium-dependen NADPH oxidase that generates superoxide, and functions as a calcium-dependent proton channel that may regulate redox-dependent processes in lymphocytes and spermatozoa. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:86 kDa

Positive control:293T and hela cells, human testis tissue, skov3, Raji and Jurkat cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE

Lysate: 40 µg

Lane 1-6: 293T cells, hela cells, human testis tissue, skov3


cells, Raji cells, Jurkat cells

Primary antibody: ml222775(NOX5 Antibody) at dilution

1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn