

Anti-HMGN5 antibody

Cat. No.	ml222786
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HMGN5 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human HMGN5
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HMGN5
Full name	high mobility group nucleosome binding domain 5
Synonyms	NSBP1; NBP-45
Swissprot	P82970

Target Background

This gene encodes a nuclear protein with similarities to the high mobility group proteins, HMG14 and HMG17, which suggests that this protein may function as a nucleosomal binding and transcriptional activating protein.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 32 kDa

Positive control: Raji, 293T, Hela and 231 cells

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-4: Raji cells, 293T cells, Hela cells, 231 cells

Primary antibody: ml222786(HMGN5 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn