

Anti-NUSAP1 antibody

Cat. No.	ml222800
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NUSAP1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human NUSAP1
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NUSAP1
Full name	nucleolar and spindle associated protein 1
Synonyms	LNP; ANKT; SAPL; BM037; NUSAP; Q0310; PRO0310p1
Swissprot	Q9BXS6

Target Background

NUSAP1 is a nucleolar-spindle-associated protein that plays a role in spindle microtubule organization. May associate with chromosomes and promote the organization of mitotic spindle microtubules around them.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 49 kDa

Positive control: HT-29, Hela and 293T cells

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-3: HT29 cells, Hela cells, 293T cells

Primary antibody: ml222800(NUSAP1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn