

Anti-RHOBTB1 antibody

Cat. No.	ml223021
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-RHOBTB1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human RHOBTB1
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	RHOBTB1
Full name	Rho-related BTB domain containing 1
Synonyms	
Swissprot	O94844

Target Background

The protein encoded by this gene belongs to the Rho family of the small GTPase superfamily. It contains a GTPase domain, a proline-rich region, a tandem of 2 BTB (broad complex, tramtrack, and bric-a-brac) domains, and a conserved C-terminal region. The protein plays a role in small GTPase-mediated signal transduction and the organization of the actin filament system.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 79 kDa

Positive control: A431 cells

Recommended dilution: 200-1000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane: A431 cells

Primary antibody: ml223021 (RHOB1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn