

Anti-S100PBP antibody

Cat. No.	ml223085
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-S100PBP rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human S100PBP
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	S100PBP
Full name	S100P binding protein
Synonyms	S100PBPR
Swissprot	Q96BU1

Target Background

This gene encodes a protein that was originally identified by its interaction with S100 calcium-binding protein P. Expression of this protein has been reported to be associated with pancreatic ductal adenocarcinoma. Alternatively spliced transcript variants have been described.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:46 kDa

Positive control:Human fetal brain tissue

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Human fetal brain tissue

Primary antibody: ml223085(S100PBP Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn