

Anti-SAA4 antibody

Cat. No.	ml223088
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SAA4 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Full length fusion protein
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SAA4
Full name	serum amyloid A4, constitutive
Synonyms	CSAA; C-SAA
Swissprot	P35542

Target Background

Major acute phase reactant. Apolipoprotein of the HDL complex.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 15 kDa

Positive control: Human plasma solution

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Human plasma solution

Primary antibody: ml223088(SAA4 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn