

订购热线: 4008-898-798

# Anti-APOBEC2 antibody

**Cat. No.** ml261294

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-APOBEC2 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human APOBEC2

ReactivityHuman, MouseContent0.2 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol APOBEC2

Full name apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 2

Synonyms ARP1; ARCD1
Swissprot Q9Y235

#### **Target Background**

Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 2 (APOBEC2) is a 224 amino acid protein that belongs to the cytidine and deoxycytidylate deaminase family. Expressed exclusively in heart and skeletal muscle, APOBEC2 is thought to be a probable C (cytidine) to U (uridine) editing enzyme. However, unlike other members of the family, such as APOBEC1, which mediates the editing of apolipoprotein (apo) B mRNA, APOBEC2 does not display any detectable apoB mRNA editing activity. Also, APOBEC2 has been shown to have low, but definite, intrinsic cytidine deaminase activity.



订购热线: 4008-898-798

## **Applications**

### Western blotting

Predicted band size:26 kDa

Positive control: Mouse muscle and skin tissue

Recommended dilution: 200-1000

Gel: 10%SDS-PAGE

Lysate: 40 µg


Lane 1-2: Mouse muscle tissue, Mouse skin tissue

Primary antibody: ml261294(APOBEC2 Antibody) at dilution

1/266.6

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



## **ELISA**

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn