

Anti-BLK antibody

Cat. No. ml221989

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-BLK rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human BLK

Reactivity Human, Mouse

Content 0.4 mg/ml

Host species Rabbit

Ig class Immunogen-specific rabbit IgG

Purification Antigen affinity purification

Target information

Symbol BLK

Full name B lymphoid tyrosine kinase



Synonyms MODY11

Swissprot P51451

Target Background

This gene encodes a nonreceptor tyrosine-kinase of the src family of proto-oncogenes that are typically involved in cell proliferation and differentiation. The protein has a role in B-cell receptor signaling and B-cell development. The protein also stimulates insulin synthesis and secretion in response to glucose and enhances the expression of several pancreatic beta-cell transcription factors.




Applications


Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 25-100





Good elisakit produceri

The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml221989(BLK Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:58 kDa


Positive control:A172 cells


Recommended dilution: 200-1000

Predicted cell location: Cytoplasm

Positive control: Human gastric cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml221989(BLK Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)



Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: A172 cells

Primary antibody: ml221989(BLK Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn