

## Anti-APOBEC3C antibody

<b>Cat. No.</b>	ml261295
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-APOBEC3C rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide of human APOBEC3C
<b>Reactivity</b>	Human
<b>Content</b>	0.3 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	APOBEC3C
<b>Full name</b>	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3C
<b>Synonyms</b>	A3C; PBI; ARP5; ARDC2; ARDC4; APOBEC1L; bK150C2.3
<b>Swissprot</b>	Q9NRW3

### Target Background

This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control.

订购热线: 4008-898-798

#### Applications

##### Western blotting

Predicted band size: 23 kDa

Positive control: PC3 cells

Recommended dilution: 500-2000

Gel: 12% SDS-PAGE

Lysate: 40 µg

Lane: PC3 cells

Primary antibody: ml261295 (APOBEC3C Antibody) at  
dilution 1/420

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



##### ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)