

Anti-SERPINB11 antibody

Cat. No. ml223136

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-SERPINB11 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human SERPINB11

Reactivity Human, Mouse

Content 0.4 mg/ml

Host species Rabbit

Ig class Immunogen-specific rabbit IgG

Purification Antigen affinity purification

Target information

Symbol SERPINB11

Full name serpin peptidase inhibitor, clade B (ovalbumin), member 11 (gene/pseudogene)



Synonyms EPIPIN; SERPIN11

Swissprot Q96P15

Target Background

SerpinB11 (serpin peptidase inhibitor, clade B (ovalbumin), member 11), also known as EPIPIN or SERPIN11, is a 392 amino acid cytoplasmic protein that belongs to the Ov-serpin subfamily and serpin family. Like other members of the serpin family, SerpinB11 has been identified as a noninhibitory intracellular protein. The gene encoding SerpinB11 maps to human chromosome 18, which houses over 300 protein-coding genes and contains nearly 76 million bases. There are a variety of diseases associated with defects in chromosome 18-localized genes, some of which include Trisomy 18 (also known as Edwards syndrome), Niemann-Pick disease, hereditary hemorrhagic telangiectasia, erythropoietic protoporphyria and follicular lymphomas.




Applications


Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 10-50





Good elisakil producere

The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml223136(SERPINB11 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:44 kDa

Positive control:Mouse lung tissue

Recommended dilution: 500-2000



Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Mouse lung tissue

Primary antibody: ml223136(SERPINB11 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn