

Anti-GIPR antibody

Cat. No.	ml261362
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-GIPR rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human GIPR
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	GIPR
Full name	gastric inhibitory polypeptide receptor
Synonyms	PGQTL2
Swissprot	P48546

Target Background

This gene encodes a G-protein coupled receptor for gastric inhibitory polypeptide (GIP), which was originally identified as an activity in gut extracts that inhibited gastric acid secretion and gastrin release, but subsequently was demonstrated to stimulate insulin release in the presence of elevated glucose. Mice lacking this gene exhibit higher blood glucose levels with impaired initial insulin response after oral glucose load. Defect in this gene thus may contribute to the pathogenesis of diabetes.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 53 kDa

Positive control: Mouse pancreas tissue

Recommended dilution: 200-1000

Gel: 10% SDS-PAGE


Lysate: 40 µg

Lane: Mouse pancreas tissue

Primary antibody: ml261362 (GIPR Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn