

Anti-TRIM35 antibody

Cat. No.	ml223349
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TRIM35 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human TRIM35
Reactivity	Human
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TRIM35
Full name	tripartite motif containing 35
Synonyms	HLS5; MAIR
Swissprot	Q9UPQ4

Target Background

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The function of this protein has not been identified.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 57 kDa

Positive control: 293T, HepG2, K562 and Hela cells

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-4: 293T cells, HepG2 cells, K562 cells, Hela cells

Primary antibody: ml223349 (TRIM35 Antibody) at dilution 1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn