

Anti-TIMM17A antibody

Cat. No.	D223354
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TIMM17A rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Full length fusion protein
Reactivity	Human, Mouse, Rat
Content	0.8 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TIMM17A
Full name	translocase of inner mitochondrial membrane 17 homolog A (yeast)
Synonyms	TIM17; TIM17A
Swissprot	Q99595

Target Background

Essential component of the TIM23 complex, a complex that mediates the translocation of transit peptide-containing proteins across the mitochondrial inner membrane.

订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm





Positive control: Human breast cancer

Recommended dilution: 25-100

Predicted cell location: Cytoplasm

Positive control: Human gastric cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml223354(TIMM17A Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: $\times 200$)

The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml223354(TIMM17A Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: $\times 200$)

Western blotting

Predicted band size: 18 kDa

Positive control: A549, K-562, Hela, NIH/3T3 and RAW264.7 cell

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE


Lysate: 40 μ g

Lane 1-5: A549, K-562, Hela, NIH/3T3 and RAW264.7 cell

Primary antibody: ml223354(TIMM17A Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn