

Anti-MYB antibody

Cat. No.	ml261474
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MYB rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human MYB
Reactivity	Human, Mouse
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MYB
Full name	v-myb avian myeloblastosis viral oncogene homolog
Synonyms	efg; Cmyb; c-myb; c-myb_CDS
Swissprot	P10242

Target Background

This gene encodes a transcription factor that is a member of the MYB family of transcription factor genes. The protein contains three domains, an N-terminal DNA-binding domain, a central transcriptional activation domain and a C-terminal domain involved in transcriptional repression. This protein plays an essential role in the regulation of hematopoiesis and may play a role in tumorigenesis. Alternative splicing results in multiple transcript variants.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 72 kDa

Positive control: HeLa and Raji cells

Recommended dilution: 200-1000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane 1-2: HeLa cells, Raji cells

Primary antibody: ml261474 (MYB Antibody) at dilution 1/633

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn