

Anti-SCARB1 antibody

Cat. No.	ml261530
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SCARB1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human SCARB1
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SCARB1
Full name	scavenger receptor class B, member 1
Synonyms	CLA1; SRB1; CLA-1; SR-BI; CD36L1; HDLQTL6
Swissprot	Q8WTV0

Target Background

The protein encoded by this gene is a plasma membrane receptor for high density lipoprotein cholesterol (HDL). The encoded protein mediates cholesterol transfer to and from HDL. In addition, this protein is a receptor for hepatitis C virus glycoprotein E2. Two transcript variants encoding different isoforms have been found for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 61 kDa

Positive control: 293T and PC3 cells

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-2: 293T cells, PC3 cells

Primary antibody: ml261530(SCARB1 Antibody) at dilution 1/260

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn