

订购热线: 4008-898-798

# Anti-TRIB2 antibody

Cat. No.	ml223414
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-TRIB2 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Full length fusion protein
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen affinity purification
Target information	oroc
Symbol	TRIB2
Full name	tribbles pseudokinase 2
	G00-
	$\checkmark$



Synonyms C5FW; TRB2; GS3955

Swissprot Q92519

#### **Target Background**

This gene encodes one of three members of the Tribbles family. The Tribbles members share a Trb domain, which is homologous to protein serine-threonine kinases, but lacks the active site lysine and probably lacks a catalytic function. The Tribbles proteins interact and modulate the activity of signal transduction pathways in a number of physiological and pathological processes.




#### Applications

#### Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human breast cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffinembedded Human breast cancer tissue using ml223414(TRIB2 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

#### Western blotting

Predicted band size:39 kDa

Positive control:231 cell and human fetal liver tissue

Good elisakii producer

Recommended dilution: 1000-5000



# Gel: 8%SDS-PAGE Lysate: 40 µg Lane 1-2: 231 cell and human fetal liver tissue Primary antibody: ml223414(TRIB2 Antibody) at dilution 1/450 Secondary antibody: Goat anti rabbit lgG at 1/8000 dilution Exposure time: 1 minute

#### ELISA

Recommended dilution: 2000-10000

## 联系电话: 4008-898-798,021-61725725

## 联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn