

订购热线: 4008-898-798

Anti-UBAP1 antibody

Cat. No. ml223463

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-UBAP1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human UBAP1

Reactivity Human, Mouse, Rat

Content0.2 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol UBAP1

Full name ubiquitin associated protein 1
Synonyms UAP; UBAP; NAG20; UBAP-1

Swissprot Q9NZ09

Target Background

This gene is a member of the UBA domain family, whose members include proteins having connections to ubiquitin and the ubiquitination pathway. The ubiquitin associated domain is thought to be a non-covalent ubiquitin binding domain consisting of a compact three helix bundle. This particular protein originates from a gene locus in a refined region on chromosome 9 undergoing loss of heterozygosity in nasopharyngeal carcinoma (NPC).



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:55 kDa

Positive control:Human placenta tissue

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane: Human placenta tissue

Primary antibody: ml223463(UBAP1 Antibody) at dilution

1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute

kDa 95 — 72 — 55 — 36 — 28 —

ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn