

Anti-TSPYL2 antibody

Cat. No.	ml261562
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TSPYL2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human TSPYL2
Reactivity	Human
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TSPYL2
Full name	TSPY-like 2
Synonyms	CDA1; CTCL; NP79; TSPX; CINAP; DENTT; SE204; HRIHFB2216
Swissprot	Q9H2G4

Target Background

This gene encodes a member of the testis-specific protein Y-encoded, TSPY-like/SET/nucleosome assembly protein-1 superfamily. The encoded protein is localized to the nucleolus where it functions in chromatin remodeling and as an inhibitor of cell-cycle progression. This protein may play a role in the suppression of tumor growth.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 79 kDa

Positive control: 231 cells

Recommended dilution: 500-2000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane: 231 cells

Primary antibody: ml261562 (TSPYL2 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn