

# Anti-YY1 antibody

**Cat. No.** ml223590

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-YY1 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human YY1

Reactivity Human, Mouse

Content 0.9 mg/ml

Host species Rabbit

Ig class Immunogen-specific rabbit IgG

Purification Antigen affinity purification

Target information

Symbol YY1

Full name YY1 transcription factor



Synonyms DELTA; NF-E1; UCRBP; INO80S; YIN-YANG-1

Swissprot P25490

# **Target Background**

YY1 is a ubiquitously distributed transcription factor belonging to the GLI-Kruppel class of zinc finger proteins. The protein is involved in repressing and activating a diverse number of promoters. YY1 may direct histone deacetylases and histone acetyltransferases to a promoter in order to activate or repress the promoter, thus implicating histone modification in the function of YY1.



# **Applications**


#### Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human brain

Recommended dilution: 25-100





Good elisakii producere

The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml223590(YY1 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: ×200)

# Western blotting

Predicted band size:45 kDa

Positive control:Hela and A549 cell lysates

Recommended dilution: 500-2000



Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-2: Hela and A549 cell lysates

Primary antibody: ml223590(YY1 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 3 seconds



### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn