

Anti-CARTPT antibody

Cat. No.	ml223668
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-CARTPT rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Full length fusion protein
Reactivity	Human, Mouse, Rat
Content	1 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen affinity purification
Target information	oroc
Symbol	CARTPT
Full name	CART prepropeptide
(-,005



Synonyms CART

Swissprot Q16568

Target Background

This gene encodes a secreted protein which is processed by prohormone/proprotein convertases to produce smaller, biologically active peptides. Expression of the transcript for this gene is regulated by certain drugs such as cocaine, and the encoded protein is thought to be involved in the regulation of appetite and stress. Mutations in this gene are associated with susceptibility to obesity.




Applications

Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human esophagus cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml223668(CARTPT Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:13 kDa

Positive control:Rat heart tissue and Mouse brain tissue

Good elisakii producer

Recommended dilution: 500-2000



Gel: 12%SDS-PAGE

Lysate: 40 µg

Lane 1-2: Rat heart tissue and Mouse brain tissue Primary antibody: ml223668(CARTPT Antibody) at dilution 1/350 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798,021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn

