

订购热线: 4008-898-798

# Anti-DSG1 antibody

**Cat. No.** ml261620

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-DSG1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human DSG1

ReactivityHuman, MouseContent0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol DSG1

Full name desmoglein 1

Synonyms DG1; DSG; CDHF4; PPKS1; SPPK1

Swissprot Q02413

#### **Target Background**

Desmosomes are cell-cell junctions between epithelial, myocardial and certain other cell types. Desmoglein 1 is a calcium-binding transmembrane glycoprotein component of desmosomes in vertebrate epithelial cells. Currently, three desmoglein subfamily members have been identified and all are members of the cadherin cell adhesion molecule superfamily. These desmoglein gene family members are located in a cluster on chromosome 18. The protein encoded by this gene has been identified as the autoantigen of the autoimmune skin blistering disease pemphigus foliaceus.



订购热线: 4008-898-798

#### **Applications**

### Western blotting

Predicted band size:114 kDa

Positive control:Human tongue cancer tissue

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE

Lysate: 40 µg

Lane: Human tongue cancer tissue

Primary antibody: ml261620(DSG1 Antibody) at dilution

1/750

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes

kDa 130— 95— 72—

## **ELISA**

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn