

Anti-CFLAR antibody

Cat. No.	ml261686
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CFLAR rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CFLAR
Reactivity	Human, Mouse
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CFLAR
Full name	CASP8 and FADD-like apoptosis regulator
Synonyms	CASH; FLIP; MRIT; CLARP; FLAME; Casper; FLAME1; c-FLIP; FLAME-1; I-FLICE; c-FLIPL; c-FLIPR; c-FLIPS; CASP8AP1
Swissprot	O15519

Target Background

The protein encoded by this gene is a regulator of apoptosis and is structurally similar to caspase-8. However, the encoded protein lacks caspase activity and appears to be itself cleaved into two peptides by caspase-8. Several transcript variants encoding different isoforms have been found for this gene, and partial evidence for several more variants exists.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 55 kDa

Positive control: LoVo and Raji cells, mouse pancreas and skeletal muscle tissue

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-4: LoVo cells, Raji cells, mouse pancreas tissue, skeletal muscle tissue

Primary antibody: ml261686(CFLAR Antibody) at dilution 1/150

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn