

Anti-HES5 antibody

Cat. No.	ml261739
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HES5 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human HES5
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HES5
Full name	hes family bHLH transcription factor 5
Synonyms	bHLHb38
Swissprot	Q5TA89

Target Background

This gene encodes a member of a family of basic helix-loop-helix transcriptional repressors. The protein product of this gene, which is activated downstream of the Notch pathway, regulates cell differentiation in multiple tissues. Disruptions in the normal expression of this gene have been associated with developmental diseases and cancer.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 18 kDa

Positive control: Mouse skin tissue, human fetal liver and fetal brain tissue

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane 1-3: Mouse skin tissue, human fetal liver tissue,
Human fetal brain tissue

Primary antibody: ml261739(HES5 Antibody) at dilution
1/100

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn