

Anti-HNRNPH1 antibody

Cat. No.	ml261750
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HNRNPH1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human HNRNPH1
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HNRNPH1
Full name	heterogeneous nuclear ribonucleoprotein H1 (H)
Synonyms	HNRPH; HNRPH1; hnRNPH
Swissprot	P31943

Target Background

This gene encodes a member of a subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA. These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some may shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNA and is very similar to the family member HNRPF. This gene may be associated with hereditary lymphedema type I.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 49 kDa

Positive control: LoVo and A375 cells, mouse brain tissue, HepG2, A431 and Hela cells

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-6: LoVo cells, A375 cells, mouse brain tissue, HepG2 cells, A431 cells, Hela cells

Primary antibody: ml261750(HNRNPH1 Antibody) at dilution 1/100

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 40 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn