

Anti-IFNA16 antibody

Cat. No.	ml261774
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-IFNA16 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human IFNA16
Reactivity	Human
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	IFNA16
Full name	interferon, alpha 16
Synonyms	IFN-alphaO
Swissprot	P05015

Target Background

Interferon, alpha 16 is produced by macrophages and has antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 22 kDa

Positive control: Riji cells

Recommended dilution: 500-2000

Gel: 8% SDS-PAGE


Lysate: 40 µg

Lane: Riji cells

Primary antibody: ml261774 (IFNA16 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn