

Anti-MAP1LC3C antibody

Cat. No.	ml261853
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MAP1LC3C rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human MAP1LC3C
Reactivity	Human
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MAP1LC3C
Full name	microtubule-associated protein 1 light chain 3 gamma
Synonyms	LC3C; ATG8J
Swissprot	Q9BXW4

Target Background

Autophagy is a highly regulated bulk degradation process that plays an important role in cellular maintenance and development. MAP1LC3C is an ortholog of the yeast autophagosome protein Atg8

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 17 kDa

Positive control: 823 and Raji cells

Recommended dilution: 200-1000

Gel: 10% SDS-PAGE


Lysate: 40 µg

Lane 1-2: 823 cells, Raji cells

Primary antibody: ml261853 (MAP1LC3C Antibody) at dilution 1/450

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn