

Anti-LLGL1 antibody

Cat. No.	ml261861
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-LLGL1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human LLGL1
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	LLGL1
Full name	lethal giant larvae homolog 1
Synonyms	DLG4; HUGL; LLGL; HUGL1; HUGL-1
Swissprot	Q15334

Target Background

This gene encodes a protein that is similar to a tumor suppressor in Drosophila. The protein is part of a cytoskeletal network and is associated with nonmuscle myosin II heavy chain and a kinase that specifically phosphorylates this protein at serine residues. The gene is located within the Smith-Magenis syndrome region on chromosome 17.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 115 kDa

Positive control: A549 cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: A549 cells

Primary antibody: ml261861(LLGL1 Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn